




OmniVac

Surface Analysis and Vacuum Technology

Inverse Photoelectron Spectroscopy System

IPES 2000



Density of unoccupied electronic states of the organic semiconductor PTCD on Ag measured with IPES 2000 detector.

The OmniVac IPES 2000 System is a highly sensitive surface analysis tool for probing the density of unoccupied electronic states between the Fermi and the Vacuum Level, which is not accessible for PES experiments.

- Solid state detector (channeltron)
- Bandpass filter for UV photon detection
- Low energy electron source with high beam current
- overall resolution 0.9 eV
- easy-to-use

www.omnivac.de

IPES 2000 System


The **OmniVac IPES 2000 System** is a highly sensitive surface analysis tool for probing the density of unoccupied electronic states between the Fermi and the vacuum level, which is not accessible for PES experiments.

The IPES 2000 system consists of an UV bandpass photon detector based on a NaCl coated channeltron in combination with a SrF₂ entrance window.

The electron source is fully UHV compatible and mounts on a DN 40 CF flange.

The thermionic emission of electrons is achieved from an indirectly heated BaO cathode.

The electron beam is characterized by a high beam current in the low energy range and a small energy distribution. Electron beam energy, beam current and spot size are independently adjustable over wide ranges.



IPES 2000 Detector

<i>Detector</i>	Solid State Detector (NaCl)
<i>Window</i>	SrF ₂ (standard), other materials possible
<i>Resolution</i>	0.9 eV (overall)
<i>Mounting flange</i>	DN 40 CF
<i>Data acquisition</i>	Preamplifier and pulse counter with USB interface
<i>Data processing</i>	Software for signal counting
<i>Remote control</i>	Electron source energy control (USB)

Low Energy Electron Source System

<i>Beam energy</i>	1 eV - 200 eV
<i>Beam current</i>	1 nA to 20 μ A
<i>Spot size</i>	0,5 mm to 5 mm at 20 mm working distance
<i>Working distance</i>	Range: 5 mm to 100 mm, 20 mm recommended
<i>Filament</i>	BaO
<i>Mounting flange</i>	DN 40 CF
<i>Power supply</i>	All necessary voltages to drive the electron source
<i>Computer/Remote control</i>	USB

- ▶ channeltron UV-Detector
- ▶ standard configuration
NaCl/SrF₂
resolution 0.9 eV
- ▶ window heating option
- ▶ other configurations possible



OmniVac
Hertelsbrunnenring 30
D-67657 Kaiserslautern / Germany

Tel +49 (0) 631 3110740
Fax +49 (0) 631 28056
sales@omnivac.de